

Photo-Therapy 4000

Bilirubintherapy in a new light. With its high reliability and low running costs, the Photo-Therapy 4000 is recommended as a highly effective piece of pediatric equipment for use in the breakdown of high concentrations of bilirubin



PT4000

PHOTO-THERAPY 4000 FEATURES

- equal light distribution
- high intensity
- quiet, no noise of a fan
- white light for diagnosis
- small dimensions

Light helps to break down the substance bilirubin. The liver of a newborn often has a low detoxication capacity. Treatment to reduce bilirubin is usually required if a newborn icterus occurs in the first day of life.

Photo-Therapy has proven to be a reliable method in the treatment of high bilirubin levels. Light transforms the bilirubin into a water soluble state by a photochemical reaction. This transformation enables precipitation by the gall bladder and kidneys.

The blue part of the spectrum, with an absorption peak of 460 nm is particularly effective for the disintegration of bilirubin. The Dräger "blue" fluorescent lights emit rays precisely in this range and are, therefore, ideally suited to the treatment. With its high reliability and low running costs, the Photo-Therapy 4000 is recommended

as an effective piece of pediatric equipment for use in the breakdown of high concentrations of bilirubin.

PRODUCT ADVANTAGES

- Small dimensions
- Low weight
- High intensity
- Quiet, no noise of a fan
- Additional white light for diagnosis. Blue and white light can be switched on separately or together
- For enhanced intensity a replacement of the 2 »white« light tubes with 2 »blue« light tubes is possible
- Stand: for the functional placement of Photo-Therapy 4000 we offer a stand on castors

TECHNICAL SPECIFICATIONS DRÄGER PHOTO-THERAPY 4000

Enim ad minim veniam

4 fluorescent tubes "blue" light	
2 fluorescent tubes "blue" light	
Power supply	120 watts

Radiant intensity in a distance of 40 cm/16 inches (US)

With 4 blue fluorescent lights	E _{Bi} = 1.6 ±0.3 mW/cm ²
With 6 blue fluorescent lights	E _{Bi} = 2.3 ±0.4 mW/cm ²

Counter for operating hours integrated

Dimensions (L x W x H)	54 x 28.5 x 13.5 cm 21 x 11 x 5 inches (US)
Weight	110–127 V: 7.2 kg/16 lbs. 230 V: 5.6 kg/12 lbs.
Recommended bulb exchange of »blue« light tubes	every 1000 hours

ORDER INFORMATION DRÄGER PHOTO-THERAPY 4000

Name	Order-No.
Photo-Therapy 4000	2M 22 310 2M 21 700 (US)
Stand on castors	2M 21 190
Fluorescent light "blue"	2M 21 010
Fluorescent light "white"	2M 21 009



www.draeger.com

**EUROPE, MIDDLE EAST, AFRICA,
LATIN AMERICA, ASIA, PACIFIC**
 Draeger Medical AG & Co. KG
 Moislinger Allee 53–55
 23558 Lübeck, Germany
 Tel +49 451 882-0
 Fax +49 451 882-2080
info@draeger.com

USA
 Draeger Medical, Inc.
 3135 Quarry Road
 Telford, PA 18969, USA
 Tel +1 215 721-5400
 Toll-free +1 866 343-2437
 Fax +1 215 723-5935
info.usa@draeger.com

CANADA
 Draeger Medical Canada, Inc.
 120 East Beaver Creek Road Suite 104
 Richmond Hill Ontario L4B 4V1, Canada
 Tel +1 905 736-3702
 Toll-free +1 866 343-2273
 Fax +1 905 763-1890
Canada.support@draeger.com

Manufacturer: Draeger Medical Systems, Inc., Telford, PA 18969, USA
 The quality management system at Draeger Medical Systems, Inc. is certified according to ISO 13485 and Annex II.3 of Directive 93/42/EEC (Medical devices).