

BabyFlow® – Nasal CPAP

Patient comfort, highly effective therapy and cost-efficiency - these are just some of the requirements fulfilled by the BabyFlow.



The Dräger nCPAP BabyFlow interface is at the cutting edge of medicine. It not only allows therapy that is beneficial to the patient but, as a non-invasive mode of ventilation, also protects against complications that may arise from intubation. The BabyFlow can also save money, as the high investment costs for additional ventilation devices do not apply when the BabyFlow is used on existing Babylog or Evita ventilators*.


PRONGS, MASKS AND BONNETS

The carefully produced silicone prongs and masks are available in different sizes. Their shapes are optimised anatomically, and ensure both safe use and an outstanding fit, despite specific facial features. The practical measuring tape makes it simple to determine the correct size.

The bonnets are all colour-coded to optimise processes and simplify handling. They are made of high-quality fabric and can be easily and reliably fastened, while still allowing fast access to the head.

GENTLE CARE

The BabyFlow not only supports the necessary flexibility in everyday clinic operations - its technical solution is also gentle on the patient. For example, the selected CPAP pressure level is as independent of the flow as possible, thus resulting in astonishingly quiet respiratory support. This means the BabyFlow is fully in accord with the philosophy of gentle care, which benefits the development of the infant.



BabyFlow® for disposable use with measuring tape

* Prerequisite here is a software that conforms with the operation of the BabyFlow, as well as the NIV/NeoFlow option for the Evita range.

Patient comfort, efficiency, and cost-effectiveness

ORDER INFORMATION

OVERVIEW OF SIZES FOR DISPOSABLE BONNETS

Size	Colour	Head circumference (max.)		Weight		Order number
		cm	inches	grams	Lbs	
S	Yellow	23	9.1	700 – 1,000	1.5 – 2.2	8418534
M	Light orange	26	10.2	1,000 – 1,300	2.2 – 2.9	8418535
L	Light salmon colour	28	11.0	1,300 – 1,600	2.9 – 3.5	8418536
XL	Lavender	30	11.9	1,600 – 2,400	3.5 – 5.3	8418537
XXL	Blue	33	13.0	2,400 – 3,000	5.3 – 6.6	8418538
XXL plus	Green	36	14.2	3,000 >	6.6 >	8418539

OVERVIEW OF SIZES FOR REUSABLE HEADBANDS

Size	Colour	Head circumference (max.)		Weight		Order number
		cm	inches	grams	Lbs	
XS	White	21	8.3	500 – 700	1.1 – 1.5	8418418
S	Yellow	23	9.1	700 – 1,000	1.5 – 2.2	8418419
M	Red	26	10.2	1,000 – 1,300	2.2 – 2.9	8418420
L	Green	28	11.0	1,300 – 1,600	2.9 – 3.5	8418421
XL	Orange	30	11.9	1,600 – 2,400	3.5 – 5.3	8418422
XXL	Blue	33	13.0	2,400 – 3,000	5.3 – 6.6	8418423
XXL plus	Black	36	14.2	3,000 >	6.6 >	8418540



Application of BabyFlow with mask



Application of BabyFlow with prong

ORDERING INFORMATION

Description	Order number
BabyFlow, disposable (pack of 20 pcs.)	8418583
BabyFlow, reusable	8418410
Prongs, XS (pack of 10 pcs.)	8418415
Prongs, S (pack of 10 pcs.)	8418605
Prongs, M (pack of 10 pcs.)	8418416
Prongs, L (pack of 10 pcs.)	8418531
Prongs, XL (pack of 10 pcs.)	8418417
NeoMask, S (pack of 10 pcs.)	8418491
NeoMask, M (pack of 10 pcs.)	8418490
NeoMask, L (pack of 10 pcs.)	8418619
Adapter for prong or NeoMask	8418411
Silicon hoses with connectors	8418414
Demo package (disposable)	8418550
Consisting of:	
BabyFlow, disposable (5 pcs.)	
1 x Prong, medium (pack of 10 pcs.)	
1 x Prong, medium plus (pack of 10 pcs.)	
1 x NeoMask, M (pack of 10 pcs.)	
1 x Cap, L (pack of 5 pcs.)	
1 x Cap, XL (pack of 5 pcs.)	
Each component of the demo package is also available separately	
Starter package (reusable)	8418435
Consisting of:	
2 x BabyFlow, reusable	
1 x Prong, small (pack of 10 pcs.)	
1 x Prong, medium (pack of 10 pcs.)	
1 x Prong, medium plus (pack of 10 pcs.)	
1 x Prong, large (pack of 10 pcs.)	
1 x Headband, S	
1 x Headband, M	
1 x Headband, L	
1 x Headband, XL	
1 x Headband, XXL	
1 x NeoMask, S (pack of 10 pcs.)	
1 x NeoMask, M (pack of 10 pcs.)	
Each component of the starter package is also available separately	




MT-6197-2006

BabyFlow, disposable




MT-4109-2006

BabyFlow, reusable



D-34337-2009

Prong family



D-34336-2009

NeoMask family



MT-046-2006

Disposable caps

HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

**REGION EUROPE CENTRAL
AND EUROPE NORTH**

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute
Technologie d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
d1mfr-contact@draeger.com

**REGION MIDDLE EAST, AFRICA,
CENTRAL AND SOUTH AMERICA**

Dräger Medical GmbH
Branch Office Dubai
Dubai Healthcare City, P.O. Box 505108
Dubai, United Arab Emirates
Tel + 971 436 24 762
Fax + 971 436 24 761
contactuae@draeger.com

REGION ASIA / PACIFIC

Dräger Medical
South East Asia Pte Ltd
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA

Dräger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free +1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

Manufacturer:

Dräger Medical GmbH
23542 Lübeck, Germany
The quality management system at
Dräger Medical GmbH is certified
according to ISO 13485, ISO 9001
and Annex II.3 of Directive 93/42/EEC
(Medical devices).