

## Disposable Breathing Circuits

Reliable, convenient and designed for safety. Because the breathing circuit is the direct interface to the patient, the integrity of your entire anesthesia or ventilation system depends on it. With Dräger's portfolio of disposable breathing circuits, you can rest assured that each product is carefully designed to work as part of a complete solution.



- Optimal fit and functionality helps to ensure the safety of the patient
- PVC-free design of gas leading components
- Single-patient use helps to prevent costly cross contamination
- Disposable convenience enhances workflow

### AN INTEGRAL PART OF THE PROCESS


Dräger has leveraged its decades of experience in anesthesia and ventilation devices, as well as reusable accessories, to develop these disposable solutions. A wide range of solutions help to ensure optimal fit and functionality. And we are committed to the continued expansion of our disposable line to optimize patient safety, no matter what the Dräger device, the treatment, or the age of the patient.

### THE SAFETY OF BIOCOMPATIBLE MATERIALS

To meet Dräger's rigorous quality standards, we carry out substantial biocompatibility testings. PVC-free materials for all gas leading components prevent complications due to the absence of DEHP. Plus, a choice of a latex or latex-free breathing bags in the anesthesia sets lets you tailor the breathing solution to the needs of the individual patient.

### HYGIENIC AND EFFICIENT

Because disposables are designed for single-patient use, they help prevent costly cross contamination. Disposables also contribute to process efficiency, so caregivers can allocate more time to other aspects of patient care.



Disposable Breathing Circuits

## ORDERING INFORMATION

Description	Order number
VentSet Basic	MP00300
VentStar® Basic (P) 150 without Luer Lock	MP00338
Anesthesia Set Basic, latex-free	MP00301
Anesthesia Set Basic, latex	MP00302
Anesthesia Set Flex, latex-free	MP00303
Anesthesia Set Flex, latex	MP00304
VentSet Flex	MP00305
VentSet Coax	MP00315
VentSet Coax with holder	MP00316
VentStar® MRI, 3 m	MP00327
VentStar® MRI (N) 300	MP00329
VentStar® Anaesthesia without Luer Lock	MP00330
VentStar® Anaesthesia (P) without Luer Lock	MP00331
VentStar® Watertrap 180 without Luer Lock	MP00334
VentStar® Watertrap 180	MP00337
VentStar® Bag Set	MP00222
VentStar® Bag Set 120	MP00380
VentStar® Bag Set 150	MP00381
VentStar® Bag Set (P) 110	MP00384



MT-5085-2006

Safe materials



MT-5006-2008

Hygienic and efficient

### HEADQUARTERS

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53–55  
23558 Lübeck, Germany

[www.draeger.com](http://www.draeger.com)

### REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH  
Moislinger Allee 53–55  
23558 Lübeck, Germany  
Tel +49 451 882 0  
Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

### REGION EUROPE SOUTH

Dräger Médical S.A.S.  
Parc de Haute  
Technologie d'Antony 2  
25, rue Georges Besse  
92182 Antony Cedex, France  
Tel +33 1 46 11 56 00  
Fax +33 1 40 96 97 20  
[dlimfr-contact@draeger.com](mailto:dlimfr-contact@draeger.com)

### REGION MIDDLE EAST, AFRICA, CENTRAL AND SOUTH AMERICA

Dräger Medical GmbH  
Branch Office Dubai  
Dubai Healthcare City, P.O. Box 505108  
Dubai, United Arab Emirates  
Tel + 971 436 24 762  
Fax + 971 436 24 761  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

### REGION ASIA / PACIFIC

Dräger Medical  
South East Asia Pte Ltd  
25 International Business Park  
#04-27/29 German Centre  
Singapore 609916, Singapore  
Tel +65 6572 4388  
Fax +65 6572 4399  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

### REGION NORTH AMERICA

Draeger Medical, Inc.  
3135 Quarry Road  
Telford, PA 18969-1042, USA  
Tel +1 215 721 5400  
Toll-free +1 800 437 2437  
Fax +1 215 723 5935  
[info.usa@draeger.com](mailto:info.usa@draeger.com)

### Manufacturer:

Dräger Medical GmbH  
23542 Lübeck, Germany  
The quality management system at  
Dräger Medical GmbH is certified  
according to ISO 13485, ISO 9001  
and Annex II.3 of Directive 93/42/EEC  
(Medical devices).